



Microprocessor UV-VIS Double Beam Spectrophotometer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 99999136670, 7835864003

USA Headquarters


30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific



BASIC MODE

To measure the Absorbance and transmittance




QUANTITATIVE

Coefficient Method
Standard Curve Up to 10 Standard sample may be used to establish a curve. Four methods for fitting a curve through the calibration points: Linear fit. Linear through zero, Square and cubic fit




DNA/PROTEIN TEST

Concentration and DNA purity are quickly and easily calculated:
Absorbance ratios: 260 nm/280 nm with optional subtracted absorbance at 320 nm. DNA concentration 62.9XA260-36.0XA280 Protein concentration-1552x4260-757.3xA280



WAVELENGTH SCAN

The wavelength scan intervals are 0.1,0.2,0.5,1,2,5 nm
High, Medium and low scan speed are available. They vary from 100 to 3600 nm/min
Wavelength are scanned from high to low so that the instrument waits at high WL. And it minimizes the degradation of UV sensitive samples.




KINETICS

Abs vs time graphs is displayed on the screen in real time wait time and measurement time up to 12 hours may be entered with time interval of 0.5,1,2,5,10,30 seconds and one min. Post-run manipulation includes re-scaling, curve tracking and selection of the part of the curve required for rate calculation. Rate is calculated using a linear regression algorithm before multiplying by the entered factor

APPLICATIONS

- Medicine/Pharmaceutical Industry
- Environment Monitoring
- Commodity Inspection
- Food Inspection
- Agricultural Chemistry
- Teaching in Colleges & Universities
- Metallurgy
- Geology
- Machine Manufacturing
- Petrochemical Industries
- Water and Waste water Labs
- Food and beverages Labs



Display (graphic LED 320x240 dots)



Soft touch keypad

Item Name	Quantity
Glass Cell	4 Nos.
Quartz cells	2 Nos.
Instruments Cover	1 No.
Software CD	1 No.
USB Cable	1 No.
Operational Manual	1 No.
Software Manual	1 No.
Software key	1 No.

Darsun Scientific

Microprocessor UV-VIS Double Beam Spectrophotometer

DSB-2704

Function

- Photometric: T%, Abs
- Quantitative: Standard Curve
- System Utility
- WL Scan (Spectrum Scan)
- Time Scan (Kinetics)
- DNA/Protein Test



Product Features

1. Double beam optical system
2. Low noise and Low stray light
3. Large LCD display, can display curve
4. High quality grating, detector and lamps
5. Data and Curve can be stored in real-time
6. Auto setting WL, auto Blank
7. Lamps can be turned on/off individually
8. Easy to change Pri-aligned lamps
9. Reinforced baseboard and bracket assures durability

*Design & specification are subject to change without any prior notice.

Model	DSB-2704
Wavelength Range	190-1100nm
Spectral Bandwidth	1nm
Optical System	Double Beam, Blazed Holographic Grating (1200 lines/mm)
Wavelength Accuracy	±0.5nm
Wavelength Repeatability	50.2nm
Wavelength Setting	Auto, Resolution 0.1 nm
Photometric Accuracy	±0.002 A (0-0.5A), ±0.003A (0.5-1A), ±0.3%T (0-100%T)
Photometric Repeatability	50.001 A (0-0.5A), 50.002A (0.5-1A), 50.2%T (0-100%T)
Stray Light	50.05%T (220/360nm)
Scan Speed	High, Medium, Low. Max.2000nm/minute
Baseline Flatness	±0.0015A
Stability	±0.001A/h (500nm,0A)
Noise	50.2%T/3min (250/500nm,0%T); 50.5%T/3min (250/500nm,100%T)
Sample Compartment	10mm Pathlength Cuvette
Detector	Silicon Photodiode
Lamps	Tungsten Lamp & Deuterium Lamp (Pre-aligned)
Display	Graphic LCD (320*240 Dots)
Keypad	30-key Alphanumeric Membrane Keypad
Printer	Mini Serial Printer; PC Printer
Power Requirements	AC90-250V, 50/60Hz
Dimension	545x468x245
Weight	18kg