

# Automatic Immunoassay Hormone Analyzer

## DS-101HA

### Product Features

1. Fully automated sample processing
2. Automatic calibration and quality control
3. Supports multiple sample types
4. Compact and space-saving unit

Automatic Immunoassay Hormone Analyzer DS-101HA can process over 200 tests per hour to ensure rapid and reliable results. It operates within an ambient temperature range of 15°C to 30°C, ensuring consistent performance. An internal thermal printer allows for immediate results printing, enhancing workflow efficiency. Our immunoassay analyzer is equipped with an intuitive LED touchscreen for easy navigation and operation.



Model No	DS-101HA
Throughput	Over 200 Tests/hour
Interface	Connected to LIS by COM/USB
Screen	LED Touch Screen
Printer	Internal Thermal Printer
Ambient Temperature	15°C to 30°C
Relative Humidity	40% to 85%
Atmospheric Pressure	86 kPa to 106 kPa
Power	Input 220V, Output DC 12V, 5A
Dimensions	230 × 175 × 305 mm
Weight	2.4 kg