

Automatic coagulometer
DS-202



Product Features

Test itmes Main test:
FIB Others:
Hep ProC ProS
Coagulation factor II,V,VII,VIII,IX,X,XI,XII
AT-III
AT RepT RVVT RVVCT RVVRT LA

Double magnetic route, stable performance , not liable to interference. 2. 7” Color touch screen 3. 100 assay parameters and 1000 test results can be saved 5. 58mm thermal printer 7. Clock CMOS chip inserted 8. Standard 232 serial interface 9. Power supply: 220V±15%, 50Hz, 80VA
Test method: double magnetic route 4. Temperature: 37°C±0.5°C 5. adding sample: by a timing linkage 6. precision: CV ≤3% 7. speed: ≤300s(depend on setting)

Model NO	DS-202
Operation	Automatic
Number of channels	4-channel
Test type	PT, APTT, TT
Configuration	Benchtop
Other characteristics	With touchscreen
Sample volume	0.2 ml (0.006763 US fl oz)
Number of sample positions	8-unit, 16 unit
Reagent volume	Min.: 0.05 ml (0 US fl oz) Max.: 0.2 ml (0.01 US fl oz)
Number of reagent positions	6 unit
Weight	5 kg (11 lb)
Length	370 mm (14.6 in)
Width	290 mm (11.4 in)
Height	120 mm (4.7 in)