

Manual colony counter

DS-A11



Product Overview

Displays average values
Bright LED lighting
Three-digit display for easy readings Adjustable pressure
Beeping noise prevents errors



Product Features

Manual colony counter DS-A11 accommodates a wide range of colony counts from 0 to 999 for accurate measurements. This colony counter fits petri dishes from 50 to 150 mm, offering versatility. It features multi-mode functionality for simultaneous counting from two petri dishes. Equipped with a pressure-sensitive touch-tip pen marker. Ensures thorough examination and accurate counting of colonies through cumulative marking. Provides nine times magnifications for easier counting.

Model No	DS-A11
Counting Range	0 to 999
Compatible Petri Dish	50 to 150 mm
Magnification	3 to 9 times
Digital Display	3 digits
Racks	2
Lighting	LED (white light)
Relative Moisture	0.8
Ambient Temperature	5 to 50°C
Protective Class	IP21
Power Consumption	20 W
Power	AC 100 to 240 V / 50 to 60 Hz
Dimensions	360 × 300 × 180 mm
Weight	4.0 Kg