

BENCHTOP PORTABLE MINI CENTRIFUGE

DS-6K/7K/10K



Product Overview

MINI portable centrifuge has three in one rotor: one rotor with three capacity 0.2ml, 0.5ml, 1.5/2ml. No need to change rotor, can centrifuge three capacities at the same. DS-6K/7K is also compatible with 6 pcr strips. We have three kinds of MINI centrifuges: DS-6K, DS-7K and DS-10K. DS-10K is high speed micro centrifuge.



Product Features

- > Compact design, light weight
- > Adjustable run time, 1s-99m59s
- > 800rpm to highest rpm is adjustable
- > Automatically stop spinning when lid is open
- > BLDC motor drive, stable and low noise
- > Can stop spin in 3 seconds, saving experiment time
- > Wide range power input: 100v to 240v
- > Three-in-one rotor: 8x1.5/2ml+6x0.5ml+2x8x0.2ml
no need to change rotor
- > Optional rotor for 6 PCR strips



Model	DS-6K/DS-7K/DS-10K
Max speed	6000/7500/10000rpm
Max RCF	2200/3450/6100Xg
Max capacity	8*1.5ml/2ml+6*0.5ml+2*8*0.2ml
Rotor	Three in one rotor; 6 PCR strip
Noise	≤52dB(A)/≤52dB(A)/≤57dB(A)
Power supply	AC 100-240V 50/60Hz
Dimension	200x170x130mm
Net weight	1.2kg