

Horizontal Electrophoresis

DS-SHE-3 Horizontal Electrophoresis Tank



Product Overview

Fast electrophoresis running speed, the whole process can be finished within 15mins.

96 sample detection in one time.

Be quicker to add samples and suitable for multi-channel pipettes.

It's more convenient to open the lid with the easy-open buttons.

Auto-switch-off when the lid is moved.



Product Features

DS-SHE-3 Technical Specifications

MODEL	DS-SHE-3
Gel size (WxL)	85x200 (mm)
No. of samples	10~364
Total buffer volume	~50 (ml)
Dimension (LxW xH)	300x150x61 (mm)
Weight	0.5 (kg)

DS-SHE-3 Packing List

Parts name	Quantity
Lid with cables	1
Body tank with electrodes	1
85x200mm Gel tray	1
Baffle	2
Combined comb frame	1
1.0mm thickness 10-wells comb	12
1.8mm thickness 18-wells comb	12