

Water Bath



Product Overview



Digital display thermostatic



High precise temperature control



One-piece stainless steel tank



High corrosion resistance

Model	DS-1-WB100 DS-1F-WB100	DS-2-WB100 DS-2F-WB100	DS-4-WB100 DS-4F-WB100	DS-6-WB100 DS-6F-WB100	DS-8F-WB100
Temperature Range	Room temp-100°C				
Power	300W	600W	800W	1200W	1500W
Accuracy			±0.1°C		
Working Chamber Size (mm)	160*160*130	300*150*140	300*300*140	450*300*130	620*300*120