

# DSCOD-100R Reactor



## Product Overview

1. Test-tube rack

2. 21 pcs digestion tubes



## Product Features

1. 21 samples at a time.

2. Settable heating temperature between 100-165 °C.

3. Settable heating time between 0-120 min.

4. Transparent plastic protection shield is equipped to protect the operator.

Model		DSCOD-100R
Measuring Unit	temperature	100-165 °C
	time	0-120 min
General	Power Supply	AC Adapter, 200-240 V AC input
	Dimensions	330×228×132mm
	Net Weight	4 kg
	Dimensions (Carton)	575X435X285 mm
	Shipping Weight	5.5 kg