

DSD-265, 100 Liter Freezer (-26~-10° C) , 165 Liters Refrigerator (2-8° C) Pharmacy Refrigerator with Freezer



Product Features

TEMPERATURE CONTROL SYSTEM

1. Microprocessor control system, refrigerator chamber and freezer chamber temperature can display at the same time;
2. Refrigerator chamber 2~ 8°C,
3. Freezer chamber -10~-26°C ;
4. According to needs, you can close the refrigerator chamber Or freezer chamber independently.

SECURITY CONTROL SYSTEM

1. Multi-malfunction alarms: high temperature alarm, low temperature alarm, sensor malfunction alarm, open the refrigerator chamber door alarm;
2. Two alarming methods: buzzer alarming, light flashing alarming;
3. Protection function: power-on delay, time between stops.
4. Refrigeration System
5. Optimized refrigeration technology makes the refrigerating capacity stronger;
6. CFC-free technology: it adopts CFC free refrigeration system and foam technology, good for environment protection;
7. Air-cooling design in the refrigerator chamber can ensure stable temperature.

HUMANIZED DESIGN

1. LCD display;
2. Refrigerator chamber equipped with stainless steel shelves and freezer chamber with basket drawer;
3. Safe door lock design prevents arbitrary open;
4. Innovative door lock design locks the refrigerator and freezer chambers separately and protects the contents inside safe;
5. Refrigerator chamber and freezer chamber can be closed separately, energy-saving and intelligent;
6. Wide voltage design suits for wide range of voltage.



Model No.	Total Volume (L)	Inner Temperature (°C)	Voltage	Packing Dimension (W×D×H) (mm)
YCD-265	265L	2~8, -10~-26	220V/50W	656*696*1739

Heat-Dissipating Method	Net Weight	Type	Over-Temperature Alarm
Direct	65	Cryobiology Cooling Device	with Alarm Function