

Rotary Rotating Mixer Shaker



Timing Function



Large handling capacity



Adjustable angle



Brushless motor with
strong power

TECHNICAL DATA

Model	DSRMO-80	DSRMO-80P
Motor	Brushless motor	Brushless motor
Control Mode	Knob	Digital display & Button
Type	Non-stop run	Non-stop run/Timing
Speed	20-80rpm	
Angle Range	40°-80°	
SPECS	Centrifuge tube : 50ml*8 or 15ml*6 or 1.5ml/2ml*24 Blood collection tube : 10ml*6 or 2ml/5ml*24	
Timing	/	1min - 23h59min
Allowed Ambient Humidity Temp.	5 ~ 40°C 80%RH	
Voltage	AC 100 - 220V, 50/60Hz	
Power	25W	
Dimension (mm)	440*255*102	