

Water Bath With Magnetic Stirrer



Product Overview



PID temperature control



Brushless Motor



Stainless steel inner tank



Strong Stirring Power



Model	DS-WB-6S	DS-WBS-6P
Speed / Increment	/	100-1200RPM / 50RPM
Temperature Range	Room Temp. +5°C -100°	Room Temp. +5°C -100°
Liner Size	200*200*130mm(depth)-2L flask or 3L beaker	
Product Size (mm)	260*265*325	260*265*325
Display Accuracy	0.1°C	0.1°C
Power	800W	
Timing Range	1min-23h59min	