

3D Oscillator **Shaker**

● ○



Timing function



Suitable for a variety of containers



Adjustable speed



3D mixing

TECHNICAL DATA

Model	DSSK-3D
Type	3D Oscillating
Diaplay	LCD
Tit Angle	10°
Max Load	0.8kg
Speed	10 - 80rpm
Timing	1min - 8h
Allowed Ambient Humidity Temp.	5~40°C 80%RH
Voltage	AC 100-240V, 50/60Hz
Power	10W
Dimension (mm)	184*184*130