

Digital Display Hotplate **Stirrer**



Strong magnetic field



Stainless steel ceramic coated
countertops

TECHNICAL DATA

Model	DSHS-8P
Control Mode	Brushless motor
Max Stirring Volume	5L(H ₂ O)
RPM	50-1500rpm
Timing Range	1min - 12h
Control Mode	Digital display & Knob & Button
Working Plate (mm)	Square Plate160*160mm
Plate Temp.	310°C
Allowable Ambient Temp	5~40°C 80%RH
Power	500W
Product Size (mm)	285*180*90