

Digital Mini Centrifuge



LCD display



Multi-rotor, more working capacity



Run smoothly, precise separation



Speed adjustable and timing

TECHNICAL DATA

Model	DSMC-5P	DSMC-7P
Control Mode	Digital display & Button	
Display	LCD display	
Max Relative Centrifugal Force	1680xg	3286xg
Max Rotation Speed	5000rpm±5%	7000rpm±5%
Rotor Capacity	0.2/0.5/1.5/2ml*12 + 0.2ml*8*4 Centrifuge tube 0.5/1.0ml*12 Micro blood tube Hematocrit rotor: 40mm*12 Capillary tube	
Timing	10sec-99min	
Power	25W	