

Low Speed Centrifuge

● ○



LCD display



Electronic door lock



Emergency switches for easier use



Program memory

TECHNICAL DATA

Model	DSL-400
Speed Range	500-4500rpm
Max Relative Centrifuge Force	2180xg
Rotor Capacity	15ml*4 Centrifuge tube 2/5/10ml*8 Blood tube
Increment	100rpm
Timing Range	30s ~ 59min50s
Noise	≤ 60dB(A)
Voltage	AC220V
Power	40W
Dimensions (mm)	265*260*190