



Can be set to run synchronously or asynchronously



Timing function



Stepless speed regulation, wide speed range



LCD display



Durable brushless motor

TECHNICAL DATA

| Model | DSMMS4P | DSMHS6P | DSMMS9P |
|---------------------------|----------------------|---------------|---------------|
| Control | Asynchronous control | | |
| Number Of Stirring Points | 4 | 6 | 9 |
| Max Capacity | 4*500ml | 6*500ml | 9*500ml |
| Working Plate Size | 4-(110*110mm) | 6-(110*110mm) | 9-(110*110mm) |
| Motor | Brushless motor | | |
| Speed | 50-1500rpm | | |
| Bar Size (mm) | Z30 | | |
| Timing | 1min-19h59min | | |
| Voltage | AC 100-220V | | |