

Low Speed Centrifuge

● ○



Program memory



Multi rotors optional

Temperature protection &
overshoot protection

Double electronic lock

TECHNICAL DATA

Model	DSL500-8	DSL500-12	DSL500-18	DSL500-24	DSL500-6
Rotor Capacity	15ml*8 Centrifuge tube	15ml*12 Centrifuge tube	15/50ml*6 Centrifuge tube	2/5/10ml*18	2/5/10ml*24
	2/5/10ml*8 Blood tube	2/5/10ml*12 Blood tube	2/5/10ml*6 Blood tube	Blood tube	Blood tube
Speed	500 - 4000rpm				
Max RCF	1967xg	2075xg	2146xg	2182xg	2164xg
Semidiameter Of Rotor	11cm	11.6 cm	12cm	12.2 cm	12.1 cm
Timing Range	30s - 59min50s				
Protection Function	Lid opening protection, cavity overtemperature protection, vibration protection				
Power	60W	60W	80W	80W	80W
Power Supply	AC 110V - 220V, 50/60Hz				