

Microplate Centrifuge

● ○



LCD display



Rotor for microplate, more working capacity



Brushless motor with strong power



Speed adjustable and timing

TECHNICAL DATA

Model	DSL-500-MP
Max Capacity	96-well PCR microplate * 4
Speed Range	500 - 2000 rpm
Max RCF	500xg
Semidiameter Of Rotor	10.2cm
Protection Function	Lid opening protection, cavity overtemperature protection, vibration protection
Timing	30s - 59min50s
Power	50W
Power Supply	AC 110V - 220V 50/60Hz