


Microplate oscillator

● ○



Timing function



Brushless motor with strong power

TECHNICAL DATA

Model	DSVM-600	DSVM-600D
Motor Type	Brushless motor	
Horizontal Gyration Amplitude	3mm	
Speed	500-1200rpm	
RPM Increment	10rpm	
Timing	1min - 8h	
Handling range	ELISA Plate / Tissue Culture Plate *4	
Allowed Ambient Humidity Temp.	5 ~ 40°C 80%RH	
Voltage	AC 110 - 220V, 50/60Hz	
Power	20W	
Dimension (mm)	355*260*95	355*260*105