

Microplate oscillator

● ○



Timing function



Brushless motor with strong power

TECHNICAL DATA

Model	DSVM-500-1P	DSVM-500-2P
Motor Type	Brushless motor	
Horizontal Gyration Amplitude	4mm	
Speed	100-1200rpm	
RPM Increment	10rpm	
Timing	1min - 8h	
Handling range	ELISA Plate / Tissue Culture Plate / 96-well plate	
Allowed Ambient Humidity Temp	5 ~ 40°C 80%RH	
Voltage	AC 100 - 240V, 50/60Hz	
Power	20W	
Dimension (mm)	155*150*90	180*150*90