

Mini Magnetic Stirrer



Strong magnetic field



Large mixing ability



Durable brushless motor



LCD display

TECHNICAL DATA

Model	DSMS3	DSMS5	DSMS5S
LCD Display	No	No	Yes
Control	Knob	Knob	Knob
Motor Type	Brushless motor		
Working Plate Size (mm)	125	135	135
Speed Range	200-1800rpm	350-1800rpm	200-2000pm
Max Stirring Volume	1L(H2O)	2L(H2O)	2L(H2O)
Stir Bar Size (mm)	30	30	30
Power	3W	5W	5W
Shell Protection Level	IP42		
Voltage	100-240V 50/60Hz		