

Mini Vortex Mixer



Stepless speed regulation
Wide speed range



High quality and colorful
silicone



Press to start / Compact
body



Brushless motor / Durable

TECHNICAL DATA

Model	DSVM-300	DSVM-300S
Control Mode	Constant speed	Knob
Motor	Brushless motor	
Amplitude	4.5mm (Circular oscillation)	
Maximum Sample Processing Capacity	50ml	
Speed Range	3000rpm	0-3000rpm
Shell Protection Level	IP42	
Power	12W	
Voltage	100-240V, 50/60Hz	
Product Size	93*93*72mm	100*93*72mm
Allowable Ambient Temperature	5-40°C,80%RH	
Gross / Net weight	0.55/0.36kg	0.56/0.37kg