

Multiple Magnetic Stirrer With Heating



LCD display



Heat more uniformly



Can be set to run synchronously or asynchronously



304 stainless steel plate, corrosion resistant



Durable brushless motor

TECHNICAL DATA

Model	DSMHS4P	DSMHS6P
Control	Asynchronous control	
Number Of Stirring Points	4	6
Max Capacity	4*500ml	6*500ml
Working Plate Size	4-(110*110mm)	6-(110*110mm)
RPM	50-1200 rpm	
Plate Material	304 Stainless steel	
Plate Temperature	120°C	
Timing	1min-19h59min	
Bar Size (mm)	Z30	Z30
Power	240W	360W
Voltage	AC 220V	