

Overhead Stirrer



Smooth start



DC permanent magnet motor



Adjustable speed



Powerful

TECHNICAL DATA

Model	DSOS-5	DSOS-10	DSOS-15S	DSOS-30S
Motor	Brushless motor			
Max Amount Of Stirring	5L	10L	10L	20L
Range Of Rotate Speed	200-1500rpm	100 - 1000rpm	100 - 2000rpm	50 - 1200rpm
Timing	/	/	1min-23h59min	1min-23h59min
Max Torque	10N.cm	20N.cm	150N.cm	30N.cm
Max Viscosity	5000mPas	10000mPas	8000mPas	16000mPas
Drill Chuck Diameter	1-10mm			
Power	20W	25W	30W	36W
Voltage	AC100-240V			