

RECIRCULATING CHILLER

DSY-50-10

MODEL NO		DSY-50-10
Capacity		10 Liter
Temperature range		-10°C-99°C
Temperature stability		± 0.5°C
Optional power		220V/50HZ
Show temperature setting way		Digital Settings for keys, LED Digital display (smallest unit: 0.1°C)
Cooling Capacity	-5°C	9348 W
	88°C	4188 W
Flow rate		30 L/mm
Practical caliber		360 mm
Lift		5-7 m
Voltage		380/50

