

Digital Display Hotplate **Stirrer**

● ○



Strong magnetic field

Stainless steel ceramic coated
countertops

TECHNICAL DATA

Model	DSHS-5C	DSHS-5P
Control Mode	Brushless motor	Brushless motor
Max Stirring Volume	2L(H ₂ O)	2L(H ₂ O)
RPM	200-1500rpm	50-1500rpm
Timing Range	/	1min-12h
Control Mode	Digital display & Knob	Digital display & Button
Working Plate (mm)	Φ135	Φ135
Plate Temp.	280°C	310°C
Liquid/Plate Temp Accuracy	±1°C/±3°C	±1°C/±3°C
Power	300W	400W
Product Size (mm)	225*155*90	225*155*90