


Dry Bath-Two Modules

● ○



LCD display



Programme memory, 9
temperature points



Overtemperature protection



Two modules, freely combined

TECHNICAL DATA

Module Options	DSDB-100-2P
Temperature Range	Room temp + 5°C ~ 100°C
Temperature Stability	±0.5°C
Temperature Homogeneity	±0.5°C
Display Accuracy	0.1°C
Heating Time	≤12min (25-100°C)
Cooling Time	(110*75*32 mm) * 2
Multi-point Operation	No
Temperature Point Can Be Set	9
Time Setting Range	5s--23h59min
Power	200W
Voltage	AC100-240V
Size	280*180*85 mm

Module Options	Numbers And Aperture Of Test Tube
DB1	96 * 0.2ml Centrifuge tube
DB2	48 * 0.5ml Centrifuge tube
DB3	35 * 1.5ml Centrifuge tube
DB4	35 * 2.0ml Centrifuge tube
DB5	24 * 5.0ml Blood tube
DB6	12 * 5 / 15ml Centrifuge tube
DB7	6 * 50ml Centrifuge tube
DB8	Flat bottom plate (ELISA plate/culture plate, etc.) (Non-replaceable modules)
DB9	Cuvettes dedicated module
DB10	Aluminum beads (Non-replaceable modules)
OEM Of Modules	Modules can be customized according to customer needs