

Mini Centrifuge



Multi-rotor, more working capacity



Low centre of gravity for stable operation



High frequency and wide voltage



Low noise

TECHNICAL DATA

Model	DSMC-4S	DSMC-7S	DSMC-4/7S
Control Mode	Constant speed		
Max Relative Centrifugal Force	1080xg	3286xg	1080/3286xg
Max Rotation Speed	4000rpm±5%	7000rpm±5%	4000/7000rpm±5%
Rotor Capacity	0.2/0.5/1.5/2ml*12 + 0.2ml*8*4 Centrifuge tube 0.5/1.0ml*12 Micro blood tube Hematocrit rotor: 40mm*12 Capillary tube		
Power	20W	25W	25W
Voltage	AC100-240V / DC24V		