


Dry Bath Shaking **Incubator**



LCD display



Multi-point operation



Overtemperature protection



Multipoint temperature operation

TECHNICAL DATA

Model	DS-DBS100C	DS-DBS100
Temperature Range	0°C ~ 100°C	15°C ~ 100°C
Temperature Stability	±0.5°C	±0.5°C
Display Accuracy	0.1°C	0.1°C
Heating Rate	5.5°C/min	5.5°C/min
Cooling Rate	2.5°C/min(100°C-Room temperature) 0.5°C/min(Below room temperature)	/
Speed / Timing	200-1500rpm / 1s-23h59min	
Mixing Diameter	3mm	
Program Memory	9 grades	
Power	200W	
Voltage	100-240V 50/60HZ	
Size	200*235*120mm	

Module Options	Numbers And Aperture Of Test Tube
DBS1	0.2ml*96 96-hole PCR microplate
DBS2	1.5ml *24 Centrifuge tube
DBS3	2ml *24 Centrifuge tube
DBS4	15ml *8 Centrifuge tube
DBS5	50ml *4 Centrifuge tube
OEM Of Modules	Modules can be customized according to customer needs