

Blood Transport Box

(Electronic Cooling) **DSD Series Features**



Product Features

- Uses semiconductor electronic cooling, semi-conductor heating, portable, easy to operate, energy-saving, and environmentally friendly.
- Outer shell is made from ABS material, with excellent strength, durability, and anti-corrosion.
- Tri-power supply mode: supports mains power (220V), car charging, and battery (standard battery 24V).
- Equipped with a high-precision DS18820 temperature sensor for accurate temperature control and real-time monitoring.
- Features digital display with precise temperature adjustment. Can set and maintain the desired temperature as needed.

Model	Capacity	Inner Size (mm)	Outer Size (mm)	Net Weight
DSD-12L	12L / 10-20 blood bags	340×190×260	420×300×290	6.6 kg
DSD-15L	15L	360×190×280	450×320×320	6.8 kg
DSD-19L	19L	380×190×250	480×320×340	8.0 kg
DSD-19L (Field Use)	19L	380×190×250	480×320×340	8.0 kg
DSD-12L	12L / 10-20 blood bags	340×190×260	420×300×290	6.6 kg