

Blood Transport Boxes

DSX Series Features



Product Features

- Made from medical-grade polyurethane integrated molding, providing excellent insulation and durability.
- Professional sealing design ensures product sealing, effectively reducing cold air loss and preventing heat ingress.
- Equipped with special insulation materials to improve temperature stability and prolong cold retention.
- Excellent structural strength, strong load-bearing capacity, and impact resistance.
- Various sizes are available to meet different usage needs.
- Can be used with ice packs for better performance.

Model	Capacity	Inner Size (mm)	Outer Size (mm)	Net Weight
DSX-6L	6L / 6 blood bags	230×195×185	330×225×280	2.8 kg
DSX-12L	12L / 12 blood bags	315×195×195	400×295×280	2.2 kg
DSX-17L	17L / 20 blood bags	340×220×240	430×285×300	3.8 kg
DSX-35L	35L / 40 blood bags	470×270×290	530×350×400	4.0 kg
DSX-53L	53L / 58 blood bags	513×315×330	610×410×450	17 kg