



Economical Freeze Dryer 12L-60



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Economical Freeze Dryer

DS12L-60A DS12L-60B DS12L-60C
DS12L-60E DS12L-60D



Model					
	DS12L-60A	DS12L-60B	DS12L-60C	DS12L-60D	DS12L-60E

condenser temperature	-60℃
Ice condenser capacity	12L/24h
Dimension (mm)	518x860x707(without chamber)
Vacuum pressure (pa)	<1Pa (standard <5Pa)
vacuum pump	2XZ-2(2L/s)
*optional vacuum pump	DRV-10(10³/h),low noisy and better performance
Refrige system	Auto Cascade
Defrost	Heating coil
Rated power	1500w
Net weight (Appro)	100kg(without vacuum pump,vacuum pump 20kg)

