



DM

Digital Mini Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Digital Mini Centrifuge

DSMC-5P/DSMC-7P



FEATURES



Multi-rotor more working capacity



Speed adjustable and timing



Run smoothly precise separation



LCD display

Model	DSMC-5P	DSMC-7P
Control Mode	Digital display & Button	
Display	LCD display	
Max Relative Centrifugal Force	1680xg	3286xg
Max Rotation Speed	5000rpm±5%	7000rpm±5%
Rotor Capacity	0.2/0.5/1.5/2ml*12 + 0.2ml*8*4 Centrifuge tube 0.5/1.0ml*12 Micro blood tube Hematocrit rotor: 40mm*12 Capillary tube	
Timing	10sec-99min	
Power	25W	