



DB

Dry Bath Shaking Incubator



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Dry Bath Shaking Incubator

DS-DB100C/DS-DB100



FEATURES



LCD display



Over temperature protection



Program memory, 9 temperature points



Two modules, freely combined

Model	DS-DBS100C	DS-DBS100
Temperature Range	0°C ~ 100°C	15°C ~ 100°C
Temperature Stability	±0.5°C	±0.5°C
Display Accuracy	0.1°C	0.1°C
Heating Rate	5.5°C/min	5.5°C/min
Cooling Rate	2.5°C/min(100°C-Room temperature) 0.5°C/min(Below room temperature)	/
Speed / Timing	200~1500rpm / 1s-23h59min	
Mixing Diameter	3mm	
Program Memory	9 grades	
Power	200W	
Voltage	100-240V 50/60HZ	
Size	200*235*120mm	

module options	Numbers and aperture of test tube
Db1	0.2ml*96 96-hole PCR microplate
Db2	1.5ml *24 centrifuge tube
Db3	2ml *24 centrifuge tube
Db4	15ml *8 centrifuge tube
Db5	50ml *4 centrifuge tube
ome of modules	Modules can be customized according to customer needs