



DB

Dry Bath-Two Modules



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Dry Bath-Two Modules

DSDB-100-2P



FEATURES



LCD display



Over temperature protection



Program memory, 9 temperature points



Two modules, freely combined

Module Options	DSDB-100-2P
Temperature Range	Room temp + 5°C ~ 100°C
Temperature Stability	±0.5°C
Temperature Homogeneity	±0.5°C
Display Accuracy	0.1°C
Heating Time	≤12min (25-100°C)
Cooling Time	(110*75*32 mm) * 2
Multi-point Operation	No
Temperature Point Can Be Set	9
Time Setting Range	5s--23h59min
Power	200W
Voltage	AC100-240V
Size	280*180*85 mm
Module Options	Numbers And Aperture Of Test Tube
DB1	96 * 0.2ml Centrifuge tube
DB2	48 * 0.5ml Centrifuge tube
DB3	35 * 1.5ml Centrifuge tube
DB4	35 * 2.0ml Centrifuge tube
DB5	24 * 5.0ml Blood tube
DB6	12 * 5 / 15ml Centrifuge tube
DB7	6 * 50ml Centrifuge tube
DB8	Flat bottom plate (ELISA plate/culture plate, etc.) (Non-replaceable modules)
DB9	Cuvettes dedicated module
DB10	Aluminum beads (Non-replaceable modules)
OEM Of Modules	Modules can be customized according to customer needs