



HS

High Speed Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

FEATURES



Temperature protection function



Multi-rotor, more working capacity



Program memory



Stable rotation

Darsun Scientific

High Speed Centrifuge

DSMC-10P DSMC-12P



Model	DSMC-10P	DSMC-12P
Control Mode	Digital display & Button	
Max Relative Centrifugal Force	6708xg	9660xg
Speed Range	1000-10000rpm±5%	1000-12000rpm±5%
Rotor Capacity	0.2/0.5/1.5/2ml*12 + 0.2ml*8*4 Centrifuge tube 0.5/1.0ml*12 Micro blood tube Hematocrit rotor: 40mm*12 Capillary tube	
Max Speed Acceleration / Reduction	≤12 s/ ≤15s	
Noise Under Max Rotation Speed	≤47db (A)	
Input Power	48W	
Shell Protection Level	AC100-240V / DC24V	
Allow Surrounding Temperature	5 - 40°C / 80%RH	