



WB

Infrared Sterilizer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Infrared Sterilizer

DS-900/DS-900S



FEATURES



One-piece stainless steel tank



High corrosion resistance



No open-flame, more safe



Easy to clean

Model	DS-900	DS-900S
Power	150W	240W
Max. Temperature Inside The Heating Hole	825°C±50°C	
Heating Temperature	Two temperatures available 400°C / 850°C	
Max. Diameter Of Sterilizer Goods (mm)	φ14	φ35
Length Of Heating Zone (mm)	150	100
Warm-up Time	10min	