



MP Microplate Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

FEATURES



**Rotor for microplate
more working
capacity**



**Brushless motor
with strong power**



**Speed adjustable
and timing**



LCD display

Darsun Scientific

Microplate
Centrifuge

DSL-500-MP



Model	DSL-500-MP
Max Capacity	96-well PCR microplate * 4
Speed Range	500 - 2000 rpm
Max RCF	500xg
Semidiameter Of Rotor	10.2cm
Protection Function	Lid opening protection, cavity over temperature protection, vibration projection
Timing	30s - 59min50s
Power	50W
Power Supply	AC 110V - 220V 50/60Hz