



VM **Mini Vortex mixer**



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

FEATURES



Stepless speed regulation
Wide speed range



Brushless motor
Durable



Press to start
Compact body



High quality and
colorful silicone

Darsun Scientific

Mini Vortex mixer

DSVM-300/DSVM-300S



Model	DSVM-300	DSVM-300S
Control Mode	Constant speed	Knob
Motor	Brushless motor	
Amplitude	4.5mm (Circular oscillation)	
Maximum Sample Processing Capacity	50ml	
Speed Range	3000rpm	0-3000rpm
Shell Protection Level	IP42	
Power	12W	
Voltage	100-240V, 50/60Hz	
Product Size	93*93*72m	100*93*72mm
Allowable Ambient Temperature	5-40°C,80%RH	
Gross / Net weight	0.55/0.36kg	0.56/0.37kg