



DB

Dry Bath Shaking Incubator



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Orbital Shaker

DSOS-40



Product Features

- 1.Speed: 40-200rpm
- 2.Digital panel with speed and time control
- 3.Working pan size 32X32cm,with max load 5KG
- 4.Continuous operation and timed operation range 1min-23h59min
- 5.DC brushless motor, safe, long life and free of maintenance. Operating silent and stable

MODEL NO	DSOS-40
Motor type	Brushless motor
Plate size	320*320
Turn	Circular
Maximum load	5 kg
Orbit diameter	20mm
Speed range	40-200 rpm
Timing range	1min-23hr59min
Voltage	100-220V,50/60HZ
Allowable	5-40°C 80% RH
Size	435*330*140mm
Net weight	8.1 kg