



## COD REACTOR



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

### Australia Headquarters

21b Lawler Drive, Oran Park, NSW,  
2570, Australia  
info@darsunscientific.com  
+61-283135514, 430583632

**Darsun**Scientific

Darsun Scientific

# REACTOR

DS-COD-100R



## GENERAL FEATURES

- 21 samples at a time.
- Settable heating temperature between 100-165 °C.
- Settable heating time between 0-120 min.
- Transparent plastic protection shield is equipped to protect the operator.

## INCLUDES

- 5 pcs turbidity bottles
- Carrying Case

Model	COD-100R	
Measuring Unit	temperature	100-165 °C
	time	0-120 min
General	Power Supply	AC Adapter, 200-240 V AC input
	Dimensions	330 × 228 × 132mm
	Net Weight	4 kg
	Dimensions (Carton)	575X435X285 mm
	Shipping Weight	5.5 kg