



RM Roller mixer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

FEATURES



Timing function



Anti - lock design, more convenient



Adjustable speed



Brushless motor with strong power

Darsun Scientific Roller mixer

DSRM-6/DSRM-6P/DSRM-8P



Model	DSRM-6	DSRM-6P	DSRM-8P
Motor Type	Brush motor	Brushless motor	Brushless motor
Display	/	LCD	LCD
Speed	20 - 80rpm		
Number Of Rollers	6	6	8
Roller Length (mm)	280		
Roller Center Distance (mm)	α34		
SPECS	φ 6 - 34mm		
Timing	/	1min - 23h59mi	1min - 23h59min
Allowed Ambient Humidity Temp.	5 ~ 40°C 80%RH		
Voltage	AC 100 - 220V, 50/60Hz		
Power	25W		
Dimension (mm)	440*255*105	440*255*105	430*320*100