



RR

Rotary Rotating Mixer shaker



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

FEATURES



Timing function



Large handling capacity



Adjustable angle



Brushless motor with strong power

Darsun Scientific

Rotary Rotating Mixer shaker

DSRMO-80/DSRMO-80P



Model	DSRMO-80	DSRMO-80P
Motor	Brushless motor	Brushless motor
Control Mode	Knob	Digital display & Button
Type	Non-stop run	Non-stop run/Timing
Speed	20-80rpm	
Angle Range	40°-80°	
SPECS	Centrifuge tube : 50ml*8 or 15ml*6 or 1.5ml/2ml*24 Blood collection tube : 10ml*6 or 2ml/5ml*24	
Timing	/	1min - 23h59min
Allowed Ambient Humidity Temp.	5 ~ 40°C 80%RH	
Voltage	AC 100 - 220V, 50/60Hz	
Power	25W	
Dimension (mm)	440*255*102	