



UC

Ultrasonic cleaning Bath



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific


Ultrasonic cleaning Bath

DS20D/DS 30D/DS 60D/DS 100D
DS 150D/DS 200D/DS 300D

Enhanced vibration, stable output
Adjustable time, stainless steel inner tank
Digital display screen, all-round cleaning
Adjustable temperature



Model	DS20D	DS 30D	DS 60D	DS 100D	DS 150D	DS 200D	DS 300D
Tank Capacity	2L	3L	6L	10L	15L	20L	28L
Operating Frequency	40 kHz	40 kHz	40 kHz	40 kHz	40 kHz	40 kHz	40 kHz
Number of Vibration Heads	1	2	3	4	6	6	8
Timing Range	0–60 min						
Settable Temperature	Room temperature – 80°C						
Ultrasonic Power	60W	120W	180W	240W	360W	360W	480W
Heating	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Degassing Function	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Drain Valve	No	No	No	Yes	Yes	Yes	Yes
Handles	No	Yes	Yes	Yes	Yes	Yes	Yes
Inner Tank Size (mm)	150×138×100	238×138×100	302×152×150	302×239×150	325×295×150	498×296×150	502×298×200
Overall Dimensions (mm)	161×149×217	251×149×217	311×161×272	308×248×315	335×305×315	510×308×343	510×308×343



DarsunScientific