



LS

Low Speed Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Low Speed Centrifuge

DSLCL-400



FEATURES



Electronic door lock



Program memory



Emergency switches
for easier use



LCD display

Model	DSLCL-400
Speed Range	500-4500rpm
Max Relative Centrifuge Force	2180xg
Rotor Capacity	15ml*4 Centrifuge tube 2/5/10ml*8 Blood tube
Increment	100rpm
Timing Range	30s ~ 59min50s
Noise	≤ 60dB(A)
Voltage	AC220V
Power	40W
Dimensions (mm)	265*260*190