



MC Mini Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

FEATURES



**Low center of gravity
for stable operation**



**High frequency
and wide voltage**



**Multi-rotor more
working capacity**



Low Noise

Darsun Scientific

Mini Centrifuge

DSMC-4S/DSMC-7S
DSMC-4/7S



Model	DSMC-4S	DSMC-7S	DSMC-4/7S
Control Mode	Constant speed		
Max Relative Centrifugal Force	1080xg	3286xg	1080/3286xg
Max Rotation Speed	4000rpm±5%	7000rpm±5%	4000/7000rpm±5%
Rotor Capacity	0.2/0.5/1.5/2ml*12 + 0.2ml*8*4 Centrifuge tube 0.5/1.0ml*12 Micro blood tube Hematocrit rotor: 40mm*12 Capillary tube		
Power	20W	25W	25W
Voltage	AC100-240V / DC24V		