



MC Mini Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

FEATURES



Stop running after
opening lid



Low noise

Darsun Scientific Mini Centrifuge

DSMC-4K/DSMC-7K/DSMC-10K/
DSMC-12K



Model	DSMC-4K	DSMC-7K	DSMC-10K	DSMC-12K
Max Speed	4000rpm	7000rpm	10000rpm	12000rpm
Max RCF	980xg	3010xg	6149xg	8855xg
Rotor Capacity	0.2ml*16 + 1.5 / 2.0ml*8		1.5 / 2.0ml*8	
Power	15W	20W	30W	40W
Dimensions (mm)	160*135*110	160*135*110	160*135*110	160*135*110
Weight (kg)	0.8	0.8	0.8	0.8