



PH

Base-Model pH Meter



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Base-Model pH Meter

DS-THE



pH Meter

Model		DS-THE
Parameters		pH/mV/°C
Range	pH	-2 ~ 18.00pH
	mv	±1999mV
	Temp	(0.0 ~ 100)°C
Resolution	PH	0.01 pH
	mv	1mV
	Temp	0.1°C
Accuracy	pH	±0.01pH±1bit pH
	mv	±0.1%FS
	Temp	±0.3°C±1bit
Temp, compensation		(0.0 ~ 100)°C, auto/manual
Calibration		2 point
Power		220V/110V
Dimension N.W		300×200×90 mm / 1.5kg
Shipping size G,W		330×260×280 mm / 2kg
Standard accessories		Electrode Holder
		pH composite electrode
		pH buffer solution
		Power adapter
		Manual
		Temperature probe