



Basic Power Supply



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Basic Power Supply

DS-300



Product Overview

DS series Basic power supply comprises 2 models, including: DS-300 and DS-600, which can meet the conventional horizontal electrophoresis, small and medium-sized vertical electrophoresis, membrane acetate electrophoresis experiment requirement, the lightweight and compact, simple and practical characteristics, make it popular in for basic teaching.

Product Features

Type of Output: Constant voltage or constant current Output range: DS-300: 5~300V, 1~300mA, maximum 90W Increment: 1V, 1mA
Timing Range: 1min ~ 9hr, 59min Display: LED
Output jacks: Two sets of output jacks Dimension (LxWxH): 280x 237x 118 (mm) Weight: 3.0 (kg)
Molding shell, touch keys, microprocessor intelligent control;
Automatic memory;
Standard timing operation;
Automatic detection of no-load, over-load

