



Basic Power Supply



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Basic Power Supply

DS-600



Product Overview

Type of Output: Constant voltage or constant current

Output range: DS600: 5~600V, 1~200mA, maximum 120W Increment: 1V, 1mA

Timing Range: 1min ~ 9hr, 59min Display: LED Output jacks: Two sets of output jacks

Dimension (LxWxH): 280x 237x 118 (mm)

Weight: 3.0 (kg)

Molding shell, touch keys, microprocessor intelligent control;

Automatic memory;

Standard timing operation;

Automatic detection of no-load, over-load

